


# OraLine®

The ONLY instant oral fluid drug screen that meets the SAMHSA cutoff for THC\*


## \*THC 4ng/mL (Marijuana)




No additional oral fluid needed for confirmation

**ORALINE® Cutoffs**

THC	4ng/mL
COC	25ng/mL
MET/AMP	50ng/mL
OPI	40ng/mL




5 green lines indicate device readiness.



After specimen collection, a red line opposite the control "C" symbol indicates the screen is valid.



Red lines of any intensity indicate a negative result (no drug detected).



Example above: The screen is negative for opiates and methamphetamine but positive for THC and cocaine.  
Total test time is 10 minutes.





Available online at www.sciencedirect.com

SCIENCE @ DIRECT®

Forensic Science International 161 (2006) 180–184

www.elsevier.com/locate/forensic

Forensic  
Science  
International

## Oral fluid testing for cannabis: On-site OraLine® IV s.a.t. device versus GC/MS

Vincent Cirimele <sup>a,\*</sup>, Marion Villain <sup>a</sup>, Patrick Mura <sup>b</sup>, Marc Bernard <sup>c</sup>, Pascal Kintz <sup>a</sup>

<sup>a</sup> Laboratoire ChemTox, 3, rue Grüniger, F-67400 Illkirch, France

<sup>b</sup> Laboratoire de biochimie, Poitiers, France

<sup>c</sup> Urgences Médico-judiciaires, Compiègne, France

Received 27 November 2005; received in revised form 11 January 2006; accepted 17 January 2006

### Abstract

Saliva or “oral fluid” has been presented as an alternative matrix to document drug use. The non-invasive collection of a saliva sample, which is relatively easy to perform and can be achieved under close supervision, is one of the most important benefits in a driving under the influence situation. Moreover, the presence of D9-tetrahydrocannabinol (THC) in oral fluid is a better indication of recent use than when 11-nor-D9-tetrahydrocannabinol-9-carboxylic acid (THC-COOH) is detected in urine, so there is a higher probability that the subject is experiencing pharmacological effects at the time of sampling. In the first part of the study, 27 drug addicts were tested for the presence of THC using the OraLine® IV s.a.t. device to establish the potential of this new on-site DOA detection technique. In parallel, oral fluid was collected with the Intercept® DOA Oral Specimen Collection device and tested for THC by gas chromatography mass spectrometry (GC/MS) after methylation for THC (limit of quantification: 1 ng/mL). The OraLine® device correctly identified nine saliva specimens positive for cannabis with THC concentrations ranging from 3 to 265 ng/mL, but remained negative in four other samples where low THC concentrations were detected by GC/MS (1–13 ng/mL). One false positive was noted. Secondly, two male subjects were screened in saliva using the OraLine® and Intercept® devices after consumption of a single cannabis cigarette containing 25 mg of THC. Saliva was first tested with the OraLine® device and then collected with the Intercept® device for GC/MS confirmation. In one subject, the OraLine® on-site test was positive for THC for 2 h following drug intake with THC concentrations decreasing from 196 to 16 ng/mL, while the test remained positive for 1.5 h for the second subject (THC concentrations ranging from 199 to 11 ng/mL). These preliminary results obtained with the OraLine® IV s.a.t. device indicate more encouraging data for the detection of THC using on-site tests than previous evaluations.

© 2006 Elsevier Ireland Ltd. All rights reserved.  
Keywords: Tetrahydrocannabinol (TCH); On-site test;  
DOA detection technique



Department of  
Industry and Resources  
Chemistry Centre (WA)

Robert Hansson  
08 9222 3041

EZI-TEST P/L  
85 Sheffield Road  
Welshpool WA 6106

Attention: Terry Snell - CEO

Dear Mr Snell

### ORALINE® SALIVA DRUG SCREEN

Subsequent to our discussions on the Laboratory's saliva drug confirmation procedures, I also wish to confirm that the Oraline® test device has met all its specifications with regard to accuracy, sensitivity and reliability.

As I have previously mentioned many similar products, and even more sophisticated instrumental on-site screening devices, have fallen well short with respect to sensitivity.

Our work to date shows that the Oraline® screen for cannabinoids clearly correlates with our GC-MS confirmation at the 4 ng/ml level. A THC saliva concentration of 4 ng/ml should adequately serve as a screening and confirmation reporting cut-off.

It should be also heartening for you to know that using our Laboratory's GC-MS procedure, the THC in saliva limit of quantitation has been determined to be 0.4 ng/ml. This represents an order of magnitude lower than an acceptable reporting cut-off. Our Laboratory has uniquely achieved such a low level of quantitation by derivatising the extracted THC with pentafluorobenzoyl chloride prior to analysis using GC-MS in *Negative Ion Chemical Ionisation mode*.

I must admit to being more than pleasantly surprised that the Oraline® product delivers what it claims given some of the poor products I have encountered.

I look forward to our association with EZI-Test as the Oraline® test device penetrates the marketplace; as indeed it must given its outstanding performance.

Yours sincerely

Robert Hansson  
Principal Chemist  
Forensic Science Laboratory



125 Hay Street, East Perth  
Western Australia 6004  
ABN 40 991 885 705

Telephone +618 9222 3177  
Facsimile +618 9325 7767  
www.ccwa.wa.gov.au

26<sup>th</sup> October 2004

SM 34-0068  
Rev. A 6/10/05